


## Sonon HD 500L Hand Held Ultrasound



- Truly portable point of care ultrasound to assist with the cannulation of a fistula or graft
- Offers wireless connectivity to provide freedom of movement and avoid cable drag
- Compatible with multiple smart devices – both iOS and android
- Ideal field of view especially for longitudinal views of the vein



- Particularly suitable for needle guidance and targeted access
- Designed to be maintenance free
- No subscription fees
- Future application upgrades provided free of charge

**Full application and clinical  
training included**

### Product Specification

- Operating Mode : BMode / Colour Flow Mode / PW/M Mode
- Frequency : 5 ~ 12MHz
- Depth : 1 ~ 6 cm
- Weight : 260g
- Dimensions : 57 x 178 x 39 33



### Expert Instruction by Julia McCarthy, RGN, BSc Hons

Our Clinical Nurse Trainer, Julia McCarthy, brings over 35 years of dialysis experience and a wealth of specialised knowledge. With qualifications including ENB 998 and ENB 136, Julia has served as a Renal Clinical Nurse Educator in the NHS, accumulating over three decades of expertise with plastic fistula cannulas and 18 years in ultrasound-guided needling.

Training follows the BRS/VASBI needling guidelines and UKKA clinical standards, supplemented by clinical references and Julia's extensive experience.



Please contact Xtra-MED for more information: **Email:** [info@xtra-med.com](mailto:info@xtra-med.com)

**Deborah Brady:** 07717 660461 **Steve Lane:** 07885 597850