

NHSSC Code	Product Code	Buttonhole Fistula Needles with Scab Picker (for arterial & venous site use)	Gauge	Needle Length mm	Tubing Length mm	Hub Type
FHD324	010DG	Buttonhole needle rotating hub 15G - 25mm needle with back eye, 300mm tubing, scab picker	15	25	300	Rotating
FHD318	002DG	Buttonhole needle rotating hub 16G - 25mm needle with back eye, 300mm tubing, scab picker	16	25	300	Rotating
FHD306	050DG	Buttonhole needle rotating hub 14G - 25mm needle with back eye, 300mm tubing, scab picker	14	25	300	Rotating
FHD313	021DG	Buttonhole needle rotating hub 15G - 25mm needle with back eye, 400mm tubing, scab picker	15	25	400	Rotating
FHD320	012DG	Buttonhole needle rotating hub 15G - 32mm needle with back eye, 300mm tubing, scab picker	15	32	300	Rotating
FHD308	074DG	Buttonhole needle rotating hub 15G - 20mm needle with back eye, 300mm tubing, scab picker	15	20	300	Rotating
FHD303	066DG	Buttonhole needle rotating hub 17G - 20mm needle with back eye, 300mm tubing, scab picker	17	20	300	Rotating

NHSCC Code	Product Code	Safety Fistula Needle (for arterial & venous site use)	Gauge	Needle Length mm	Tubing Length mm	Hub Type
FHD299	010SG	Safety needle rotating hub 15G - 25mm needle with back eye, protector sheath, 300mm tubing	15	25	300	Rotating
FHD296	002SG	Safety needle rotating hub 16G - 25mm needle with back eye, protector sheath, 300mm tubing	16	25	300	Rotating
FHD300	074SG	Safety needle rotating hub 15G - 20mm needle with back eye, protector sheath, 300mm tubing	15	20	300	Rotating
FHD307	070SG	Safety needle rotating hub 16G - 20mm needle with back eye, protector sheath, 300mm tubing	16	20	300	Rotating
FHD319	066SG	Safety needle rotating hub 17G - 20mm needle with back eye, protector sheath, 300mm tubing	17	20	300	Rotating
FHD342	021SG	Safety needle rotating hub 15G - 25mm needle with back eye, protector sheath, 400mm tubing	15	25	400	Rotating
FHD317	012SG	Safety needle rotating hub 15G - 32mm needle with back eye, protector sheath, 300mm tubing	15	32	300	Rotating

NHSSC Code	Product Code	Plastic Fistula Cannula Needles (soft needles for arterial & venous site use)	Gauge	Needle Length mm
FHD331	SP802-14 (25)	Supercath NEO soft safety fistula cannula 14G - 25mm cannula, protector sheath	14	25
FHD335	SP802-15 (25)	Supercath NEO soft safety fistula cannula 15G - 25mm cannula, protector sheath	15	25
FHD336	SP802-16 (25)	Supercath NEO soft safety fistula cannula 16G - 25mm cannula, protector sheath	16	25
FHD85012	SP402-14 (25)	Supercath PRO soft safety fistula cannula 14G - 25mm cannula, protector sheath, non-return valve	14	25
FHD85013	SP402-14 (33)	Supercath PRO soft safety fistula cannula 14G - 33mm cannula, protector sheath, non-return valve	14	33
FHD85016	SP402-15 (25)	Supercath PRO soft safety fistula cannula 15G - 25mm cannula, protector sheath, non-return valve	15	25
FHD85011	SP402-15 (33)	Supercath PRO soft safety fistula cannula 15G - 33mm cannula, protector sheath, non-return valve	15	33
FHD85014	SP402-16 (25)	Supercath PRO soft safety fistula cannula 16G - 25mm cannula, protector sheath, non-return valve	16	25
FHD85015	SP402-16 (33)	Supercath PRO soft safety fistula cannula 16G - 33mm cannula, protector sheath, non-return valve	16	33
FHD85009	SP402-17 (25)	Supercath PRO soft safety fistula cannula 17G - 25mm cannula, protector sheath, non-return valve	17	25
FHD85010	SP402-17 (33)	Supercath PRO soft safety fistula cannula 17G - 33mm cannula, protector sheath, non-return valve	17	33