

Product Code	Plastic Fistula Cannula Needles (soft needles for arterial & venous site use)	Gauge	Needle Length mm
SP802-14 (25)	Supercath NEO soft safety fistula cannula 14G - 25mm cannula, protector sheath	14	25
SP802-15 (25)	Supercath NEO soft safety fistula cannula 15G - 25mm cannula, protector sheath	15	25
SP802-16 (25)	Supercath NEO soft safety fistula cannula 16G - 25mm cannula, protector sheath	16	25
SP402-14 (25)	Supercath PRO soft safety fistula cannula 14G -25mm cannula, protector sheath, non-return valve	14	25
SP402-14 (33)	Supercath PRO soft safety fistula cannula 14G - 33mm cannula, protector sheath, non-return valve	14	33
SP402-15 (25)	Supercath PRO soft safety fistula cannula 15G - 25mm cannula, protector sheath, non-return valve	15	25
SP402-15 (33)	Supercath PRO soft safety fistula cannula 15G - 33mm cannula, protector sheath, non-return valve	15	33
SP402-16 (25)	Supercath PRO soft safety fistula cannula 16G - 25mm cannula, protector sheath, non-return valve	16	25
SP402-16 (33)	Supercath PRO soft safety fistula cannula 16G - 33mm cannula, protector sheath, non-return valve	16	33
SP402-17 (25)	Supercath PRO soft safety fistula cannula 17G - 25mm cannula, protector sheath, non-return valve	17	25
SP402-17 (33)	Supercath PRO soft safety fistula cannula 17G - 33mm cannula, protector sheath, non-return valve	17	33