



Blunt Buttonhole Fistula Needles With Scab Picker

Product Code	Description	Needle Gauge mm	Hub type	Back eye	Needle Length mm	Tubing Length mm	Quantity Per Case
AVF-D-BF2530BGDS-14	14G blunt buttonhole fistula needle with back eye, fixed hub, needle length 25mm, tubing length 30cm	14	Fixed	Yes	25	300	500
AVF-D-BR2530BGDS-15	15G blunt buttonhole fistula needle with back eye, rotating hub, needle length 25mm, tubing length 30cm	15	Rotating	Yes	25	300	500
AVF-D-BR2540BGDS-15	15G blunt buttonhole fistula needle with back eye, rotating hub, needle length 25mm, tubing length 40cm	15	Rotating	Yes	25	400	500
AVF-D-BR2530BGDS-16	16G blunt buttonhole fistula needle with back eye, rotating hub, needle length 25mm, tubing length 30cm	16	Rotating	Yes	25	300	500
AVF-D BR2530BGDS-17	17G blunt buttonhole fistula needle with back eye, rotating hub, needle length 25mm, tubing length 30cm	17	Rotating	Yes	25	300	500
AVF-D-BR1530BGDS-15	15G blunt buttonhole fistula needle with back eye, rotating hub, needle length 15mm, tubing length 30cm	15	Rotating	Yes	15	300	500
AVF-D-BR1530BGDS-17	17G blunt buttonhole fistula needle with back eye, rotating hub, needle length 15mm, tubing length 30cm	17	Rotating	Yes	15	300	500
AVF-D-BR3230BGDS-15	15G blunt buttonhole fistula needle with back eye, rotating hub, needle length 32mm, tubing length 30cm	15	Rotating	Yes	32	300	500

Sharp Fistula Needle With Safety Protector

Product Code	Description	Needle Gauge mm	Hub type	Back eye	Needle Length mm	Tubing Length mm	Quantity Per Case
AVF-BF2530BGP-14	14G sharp fistula needles with safety protector with backeye, fixed hub, needle length 25mm, tubing length 30cm	14	Fixed	Yes	25	300	500
AVF-BF2530BGP-15	15G sharp fistula needles with safety protector with backeye, fixed hub, needle length 25mm, tubing length 30cm	15	Fixed	Yes	25	300	500
AVF-BF2530BGP-16	16G sharp fistula needles with safety protector with backeye, fixed hub, needle length 25mm, tubing length 30cm	16	Fixed	Yes	25	300	500
AVF-BF2530BGP-17	17G sharp fistula needles with safety protector with backeye, fixed hub, needle length 25mm, tubing length 30cm	17	Fixed	Yes	25	300	500
AVF-BR2530BGP-15	15G sharp fistula needles with safety protector with backeye, rotating hub, needle length 25mm, tubing length 30cm	15	Rotating	Yes	25	300	500
AVF-BR2530BGP-16	16G sharp fistula needles with safety protector with backeye, rotating hub, needle length 25mm, tubing length 30cm	16	Rotating	Yes	25	300	500
AVF-BR2530BGP-17	17G sharp fistula needles with safety protector with backeye, rotating hub, needle length 25mm, tubing length 30cm	17	Rotating	Yes	25	300	500
AVF-BF2540BGP-15	15G sharp fistula needles with safety protector with backeye, fixed hub, needle length 25mm, tubing length 40cm	15	Fixed	Yes	25	400	500
AVF-BF3230BGP-15	15G sharp fistula needles with safety protector with backeye, fixed hub, needle length 32mm, tubing length 30cm	15	Fixed	Yes	32	300	500
AVF-BF3230BGP-16	16G sharp fistula needles with safety protector with backeye, fixed hub, needle length 32mm, tubing length 30cm	16	Fixed	Yes	32	300	500
AVF-BF3230BGP-17	17G sharp fistula needles with safety protector with backeye, fixed hub, needle length 32mm, tubing length 30cm	17	Fixed	Yes	32	300	500
AVF-BF3240BGP-15	15G sharp fistula needles with safety protector with backeye, fixed hub, needle length 32mm, tubing length 40cm	15	Fixed	Yes	32	400	500
AVF-BR3230BGP-15	15G sharp fistula needles with safety protector with backeye, rotating hub, needle length 32mm, tubing length 30cm	15	Rotating	Yes	32	300	500
AVF-BF1530BGP-15	15G sharp fistula needles with safety protector with backeye, fixed hub, needle length 15mm, tubing length 30cm	15	Fixed	Yes	15	300	500
AVF-BF1530BGP-17	17G sharp fistula needles with safety protector with backeye, fixed hub, needle length 15mm, tubing length 30cm	17	Fixed	Yes	15	300	500



Thornecroft Farm - Kettleshulme, Cheshire SK23 7RG

Full Name of Company

XtraMed, Limited (Reg. No. 5968756)

Address including Postcode

**Thornecroft Farm
Kettleshulme, Cheshire SK23 7RG**

Sales telephone number

01663 733 522 / 0779 1699 742

Telephone and fax Number for placing orders

01663 733 522

email address for placing orders

order@xtra-med.com / click99@btconnect.com

Telephone and Fax number for Remittances

01663 733 522

Name of bank

Santander UK

Sort code

09 01 27

Account number

75481418

IBAN

GB78 ABBY 0901 2775 481418

VAT number

928 6213 13

Which method is preferred for ordering (fax or email) – **email is preferred for Orders and Remittances**