

Blunt Buttonhole Fistula Needles With Scab Picker

Product Code	Needle Gauge mm	Hub type	Back eye	Suitable for both Arterial and venous sites	Needle Length mm	Tubing Length mm	Quantity Per Case
AVF-D-BR2530BGDS-14	14	Rotatable	Yes	Yes	25	300	500
AVF-D-BR2530BGDS-15	15	Rotatable	Yes	Yes	25	300	500
AVF-D-BR2540BGDS-15	15	Rotatable	Yes	Yes	25	400	500
AVF-D-BR2530BGDS-16	16	Rotatable	Yes	Yes	25	300	500
AVF-D BR2530BGDS-17	17	Rotatable	Yes	Yes	25	300	500

Sharp Safety Fistula Needles With Protector

Product Code	Needle Gauge mm	Hub type	Back eye	Suitable for both Arterial and venous sites	Needle Length mm	Tubing Length mm	Quantity Per Case
AVF-BF2530BGP-14	14	Fixed	Yes	Yes	25	300	500
AVF-BF2530BGP-15	15	Fixed	Yes	Yes	25	300	500
AVF-BF2530BGP-16	16	Fixed	Yes	Yes	25	300	500
AVF-BF2530BGP-17	17	Fixed	Yes	Yes	25	300	500
AVF-BF2540BGP-15	15	Fixed	Yes	Yes	25	400	500
AVF-BR2530BGP-15	15	Rotatable	Yes	Yes	25	300	500
AVF-BR2530BGP-16	16	Rotatable	Yes	Yes	25	300	500
AVF-BR2530BGP-17	17	Rotatable	Yes	Yes	25	300	500
AVF-BF1530BGP-15	15	Fixed	Yes	Yes	15	300	500
AVF-BF1530BGP-17	17	Fixed	Yes	Yes	15	300	500
AVF-BF1515BGP-17	17	Fixed	Yes	Yes	15	150	500
AVF-BR2030BGP-14	14	Rotatable	Yes	Yes	20	300	500
AVF-BR2030BGP-15	15	Rotatable	Yes	Yes	20	300	500
AVF-BR2030BGP-16	16	Rotatable	Yes	Yes	20	300	500
AVF-BF3230BGP-15	15	Fixed	Yes	Yes	32	300	500
AVF-BF3230BGP-16	16	Fixed	Yes	Yes	32	300	500
AVF-BF3230BGP-17	17	Fixed	Yes	Yes	32	300	500
AVF-BF3240BGP-15	15	Fixed	Yes	Yes	32	400	500
AVF-BR3230BGP-15	15	Rotatable	Yes	Yes	32	300	500