

Product Code	Buttonhole Fistula Needles With Scab Picker (for arterial & venous site use)	Gauge	Needle Length mm	Tubing Length mm	Hub Rotating / Fixed
010DG	Buttonhole needle rotating hub 15G - 25mm needle with back eye, 300mm tubing, scab picker	15	25	300	Rotating
002DG	Buttonhole needle rotating hub 16G - 25mm needle with back eye, 300mm tubing, scab picker	16	25	300	Rotating
050DG	Buttonhole needle rotating hub 14G - 25mm needle with back eye, 300mm tubing, scab picker	14	25	300	Rotating
021DG	Buttonhole needle rotating hub 15G - 25mm needle with back eye, 400mm tubing, scab picker	15	25	400	Rotating
012DG	Buttonhole needle rotating hub 15G - 32mm needle with back eye, 300mm tubing, scab picker	15	32	300	Rotating
068DG	Buttonhole needle rotating hub 17G - 15mm needle with back eye, 300mm tubing, scab picker	17	15	300	Rotating
076DG	Buttonhole needle rotating hub 15G - 15mm needle with back eye, 300mm tubing, scab picker	15	15	300	Rotating
006DG	Buttonhole needle rotating hub 17G - 25mm needle with back eye, 300mm tubing, scab picker	17	25	300	Rotating
074DG	Buttonhole needle rotating hub 15G - 20mm needle with back eye, 300mm tubing, scab picker	15	20	300	Rotating
066DG	Buttonhole needle rotating hub 17G - 20mm needle with back eye, 300mm tubing, scab picker	17	20	300	Rotating
033DG	Buttonhole needle fixed hub 15G - 25mm needle with back eye, 400mm tubing, scab picker	15	25	400	Fixed
009DG	Buttonhole needle fixed hub 15G - 25mm needle with back eye, 300mm tubing, scab picker	15	25	300	Fixed
Product Code	Safety Fistula Needle (for arterial & venous site use)	Gauge	Needle Length mm	Tubing Length mm	Hub Rotating / Fixed
009SG	Safety needle fixed hub 15G - 25mm needle with back eye, protector sheath, 300mm tubing	15	25	300	Fixed
001SG	Safety needle fixed hub 16G - 25mm needle with back eye, protector sheath, 300mm tubing	16	25	300	Fixed
005SG	Safety needle fixed hub 17G - 25mm needle with back eye, protector sheath, 300mm tubing	17	25	300	Fixed
011SG	Safety needle fixed hub 15G - 32mm needle with back eye, protector sheath, 300mm tubing	15	32	300	Fixed
023SG	Safety needle fixed hub 15G - 32mm needle with back eye, protector sheath, 400mm tubing	15	32	400	Fixed
067SG	Safety needle fixed hub 17G - 15mm needle with back eye, protector sheath, 300mm tubing	17	15	300	Fixed
069SG	Safety needle fixed hub 16G - 20mm needle with back eye, protector sheath, 300mm tubing	16	20	300	Fixed
073SG	Safety needle fixed hub 15G - 20mm needle with back eye, protector sheath, 300mm tubing	15	20	300	Fixed
010SG	Safety needle rotating hub 15G - 25mm needle with back eye, protector sheath, 300mm tubing	15	25	300	Rotating
002SG	Safety needle rotating hub 16G - 25mm needle with back eye, protector sheath, 300mm tubing	16	25	300	Rotating
006SG	Safety needle rotating hub 17G - 25mm needle with back eye, protector sheath, 300mm tubing	17	25	300	Rotating
074SG	Safety needle rotating hub 15G - 20mm needle with back eye, protector sheath, 300mm tubing	15	20	300	Rotating
070SG	Safety needle rotating hub 16G - 20mm needle with back eye, protector sheath, 300mm tubing	16	20	300	Rotating
066SG	Safety needle rotating hub 17G - 20mm needle with back eye, protector sheath, 300mm tubing	17	20	300	Rotating
050SG	Safety needle rotating hub 14G - 25mm needle with back eye, protector sheath, 300mm tubing	14	25	300	Rotating
021SG	Safety needle rotating hub 15G - 25mm needle with back eye, protector sheath, 400mm tubing	15	25	400	Rotating
012SG	Safety needle rotating hub 15G - 32mm needle with back eye, protector sheath, 300mm tubing	15	32	300	Rotating
Product Code	Plastic Fistula Cannula Needles (Soft needles for arterial & venous site use)	Gauge	Needle Length mm	Tubing Length mm	Hub Rotating / Fixed
SP802-14 (25)	Supercath Neo soft safety fistula cannula 14G, 25mm cannula, protector sheath	14	25	0	Rotating
SP802-14 (33)	Supercath Neo soft safety fistula cannula 14G, 33mm cannula, protector sheath	14	33	0	Rotating
SP802-15 (25)	Supercath Neo soft safety fistula cannula 15G, 25mm cannula, protector sheath	15	25	0	Rotating
SP802-15 (33)	Supercath Neo soft safety fistula cannula 15G, 33mm cannula, protector sheath	15	33	0	Rotating
SP802-16 (25)	Supercath Neo soft safety fistula cannula 16G, 25mm cannula, protector sheath	16	25	0	Rotating